



height adjustable
LABORATORY SEATING

WWW.BONLABIN.COM




BIMOS is Europe’s leading manufacturer of the best industrial and laboratory chairs.

The specially developed hygienic design, targeted support for positions typically found in laboratory work and the compact construction of BIMOS laboratory chairs guarantee optimal conditions in a state-of-the-art working environment.



Labsit



Models and materials

In the laboratory, the perfect cleaning process starts with the material. Labsit’s back rest is made from faux leather with 3D flex function. Depending on requirements, Labsit can be supplied with integral foam or faux leather padding for laboratories and material or Supertec padding for other uses. The characteristic seat shell comes in five attractive colours. The back rest is available in black only. The cross-shaped base is available in black plastic or polished aluminium. All of the plastic and PU parts of the stool are black. The 5-star base is also available in black plastic or polished aluminium. The optional Lab-Clip makes it easy to attach individual labels to the chair. The ESD-compliant Labsit model comprises a cross-shaped base made from polished aluminium, a black seat shell and four padding material options.

Stool functions



Functions



Labsit stool with castors

Seat height adjustable from 450mm to 610mm.

Design	Order no.
Integral foam black	9127-2000
Integral foam ESD black	9127E-2000

Options



ESD
(Seat shell =
Black)



Polished aluminium
5-star base
(ESD standard)

Accessories (for precise details, see page 149)



Multi-function
armrest 9129



Labsit 2
with foot ring



Lab-Clip
9128

Seat padding and shell

	Magic syn- thetic leather	Integral foam	Fabric Lucia	Supertec
Seat padding				
Order no.	M001 / ESD M001	2000 / ESD 2000	5800 / ESD fabric: Duotec 9801	SP01 / ESD CP01
Seat shell				
Order no.	3277 Blue	3285 Anthra- cite	3403 White	3279 Orange
				3280 Green
				3218 ESD Black



Labsit 2
with castors

Seat height adjustable from 450mm to 610mm.

Design	Order no.
Magic synthetic leather black	9123-M001-colour seat shell
Integral foam black	9123-2000-colour seat shell
Fabric Lucia black	9123-5800-colour seat shell
Supertec black	9123-SP01-colour seat shell
Magic synthetic leather ESD black	9123E-M001-3218
Integral foam ESD black	9123E-2000-3218
Fabric Duotec ESD black	9123E-9801-3218
Supertec ESD black	9123E-CP01-3218



Labsit 3
with glides and foot ring

Seat height adjustable from 520mm to 770mm.

Design	Order no.
Magic synthetic leather black	9121-M001-colour seat shell
Integral foam black	9121-2000-colour seat shell
Fabric Lucia black	9121-5800-colour seat shell
Supertec black	9121-SP01-colour seat shell
Magic synthetic leather ESD black	9121E-M001-3218
Integral foam ESD black	9121E-2000-3218
Fabric Duotec ESD black	9121E-9801-3218
Supertec ESD black	9121E-CP01-3218



Labsit 4
with Stop & Go castors and foot ring

Seat height adjustable from 560mm to 810mm.

Design	Order no.
Magic synthetic leather black	9125-M001-colour seat shell
Integral foam black	9125-2000-colour seat shell
Fabric Lucia black	9125-5800-colour seat shell
Supertec black	9125-SP01-colour seat shell
Magic synthetic leather ESD black	9125E-M001-3218
Integral foam ESD black	9125E-2000-3218
Fabric Duotec ESD black	9125E-9801-3218
Supertec ESD schwarz	9125E-CP01-3218



BLC01CW
PU SEATED LAB CHAIR WITH CASTORS

BLC01GB
PU SEATED LAB CHAIR WITH GLIDER BASE

- Height Adjustable
- Back rest tilt adjustment
- Black knob for height adjustment
- All PU foam on the back and seat
- Plastic moulded foot ring
- Chrome gas-lift
- PU Cushion and back Color: Black
- Seat Height Range: 580~835mm
- Seat Width: 450mm
- Seat Depth: 430mm
- Comfortable Lumbar Support
- PU Foam thickness: 35mm
- Base: 640mm 5-star PU base / Chrome base



BLS01CW
PU SEATED LAB STOOL WITH CASTORS

BLS01GB
PU SEATED LAB STOOL WITH GLIDER BASE

- Height Adjustable
- Black knob for height adjustment
- All PU foam for the seat
- Plastic moulded foot ring
- Chrome gas-lift
- Seat Height Range: 580~835mm
- Seat Width: 300mm
- Seat Depth: 300mm
- Comfortable Lumbar Support
- PU Foam thickness: 35mm
- Base: 640mm 5-star PU base / Chrome base





BLS02CW
WOOD SEATED LAB STOOL WITH CASTORS

BLS02GB
WOOD SEATED LAB STOOL WITH GLIDER BASE

- Height Adjustable
- Black knob for height adjustment
- Laminated / veneered / painted plywood for the seat
- Plastic moulded foot ring
- Chrome gas-lift
- Seat Height Range: 580~835mm
- Seat dia: 300mm
- Seat thickness: 24mm
- Base: 640mm 5-star PU base / Chrome base





BLS03EC
WOOD SEATED LAB STOOL FOR STUDENTS

Custom made
Height Adjustable - threaded rod
Heavy duty threaded rod for height adjustment
Laminated / veneered / painted plywood for the seat
MS Tube epoxy powder coated frame & foot ring
Seat Height Range: 580~835mm
Seat dia: 300mm
Seat thickness: 24mm



PU PAINTED FINISH - ANY RAL COULRS



FOAM CUSHIONED SEAT - SELECTED COULRS



VENEERED WOOD SEAT - SELECTED WOOD FINISH





BLC04CW
BLOW MOLDED SEATED LAB CHAIRL WITH CASTORS



BLC04GB
BLOW MOLDED SEATED LAB STOOL WITH GLIDER BASE

Height Adjustable
Black knob for height adjustment
Single piece Blow Molded Shell for seat and back rest
Plastic moulded foot ring
Chrome gas-lift
Seat Height Range: 580~835mm
Seat dia: 300mm
Seat thickness: 24mm
Base: 640mm 5-star PU base / Chrome base
Single piece Shell 911 color - Blue or Orange

